

1800 Jackson street Suite #214

Golden, co 80401

tel:1-800-537-2845

DISFIGURED EYES ORDER FORM (Hand Painted)

ACCOUNT NAME _____ ACCOUNT NUMBER _____

ADDRESS: _____ Address change _____

City: _____ State: _____ Zip: _____

Telephone: _____ Fax: _____

Patient's Name: _____ Check if previous Patient _____

WE OFFER 2 OPTIONS FOR HAND PAINTED IRISES

#1 Take a picture of both eyes together, close up in normal lighting, from nose to forehead only. We cannot use Polaroid pictures. Send the pictures to us via mail, or email. We will send you 2 Hand Painted trial lenses as close to the picture as possible. After trying on your patient, return the lenses and a picture with the closest trial lens (on the patient). If the trial lens needs modification we will continue to adjust the lens until your patient approves it. Then we will paint the patient's new lens to match the approved sample. Your patient's approved trial lens will be filed and identified for future duplications. The completed process normally takes 5 -7 weeks, however it could take longer, depending on the number of modifications. Trial lenses must be returned. !!!TRIAL LENSES ARE FOR ONLY COLOR AND PATTERN VERIFICATION - THEY MIGHT NOT BE THE RIGHT FIT FOR YOUR PATIENT !!!!!!!

#2 If your patient prefers, he /she may come to our facility in Golden, Colorado and should plan to spend 3 working days while we produce a satisfactory control trial lens. This must be done by appointment only. OFFICE VISIT FEE MUST BE PAID BY PATIENT AT THE TIME OF THE APPOINTMENT!!!!!!

Contact Lens Specifications — Please Specify the Following Information:

Orders will not be processed without a good picture and completed order form.

Your Patient has: Vision _____ No Vision _____ Light Sensitivity _____

Patient's disfigured eye is _____ Right _____ Left Description: White / Black / Other

Lens Brand and Specifications: _____

Check Pupil or Clear Center Size using attached Pupil Diameter Gauge:

2.5mm _____ 3.0mm _____ 3.5mm _____ 4.0mm _____ 4.5mm _____ 5.0mm _____ 5.5mm _____

Check Iris Size using attached Iris Diameter Gauge:

10.0mm _____ 10.5mm _____ 11.0mm _____ 11.5mm _____ 12.0mm _____ 12.5mm _____ 13.0mm

Does the Hand Painted Iris Lens Need :

White Sclera ___ Clear Sclera ___ Clear Center ___ Black Pupil ___ Black Backing ___ DECENTRATION \$50.00 EXTRA

NOTE: If design and color needs to be decentered, in order to make the eye look centered please call for instructions. Most 55% water daily wear lenses are suitable (methafilcon , ocufilecon or bufilecon material) for painting the final lens.

Warranty

COLOR ON LENS IS WARRANTED for 6 months against color fading. PLEASE CALL FOR SPECIFIC RULES !!!!!

* Lenses over 6 months old must be in good condition in order to re-color *

Order Authorized By: _____

Remarks: _____

6 month Warranty from the date of ship

PLEASE CALL FOR PRICING!!!!!!

PRICING IS FOR ONLY TINTING AND HAND PAINTING - PRICING DOES NOT INCLUDE THE LENS

*FINAL HAND PAINTED LENSES SHOULD BE DISPENSED RIGHT AWAY , SO WE CAN MAKE MODIFICATION WITHIN 30 DAYS GRACE PERIOD !!!!!!! **

** ALL ORDERS ARE CUSTOM TO MAKE - NO RETURN NOR EXCHANGE **